

## **Biomedical Instrumentation**

Prof. Dr. Nizamettin AYDIN

**Biopotential Electrodes** 





$$\begin{split} V_{\rm p} &= E^{\,0} + V_{\rm r} + V_{\rm c} + V_{\rm a} \\ E &= -\frac{RT}{nF} \ln \frac{a_1}{a_2} \\ E &= -\frac{RT}{nF} \ln (a_{c^{n+}}) \\ \alpha A + \beta B &\rightleftharpoons \gamma C + \delta D + n {\rm e}^- \\ E &= E^0 = \frac{\mu_+ - \mu_-}{\mu_+ + \mu_-} \frac{RT}{nF} \ln \frac{a'}{a''} \\ E_j &= \frac{\mu_+ - \mu_-}{\mu_+ + \mu_-} \frac{RT}{nF} \ln \frac{a'}{a''} \end{split}$$

$$Ag \rightleftharpoons Ag^{+} + e^{-}$$

$$Ag^{+} + Cl^{-} \rightleftharpoons AgCl \downarrow$$

$$a_{Ag^{+}} \times a_{Cl^{-}} = K_{s}$$

$$E = E_{Ag}^{0} + \frac{RT}{nF} \ln a_{Ag^{+}}$$

$$E = E_{Ag}^{0} + \frac{RT}{nF} \ln \frac{K_{s}}{a_{Cl^{-}}}$$

$$E = E_{Ag}^{0} + \frac{RT}{nF} \ln K_{s} - \frac{RT}{nF} \ln a_{Cl^{-}}$$

$$I = 100 \,\mathrm{mA} \, e^{-t/10 \,\mathrm{s}}$$

$$q = \int_0^\infty i \, dt = 100 \,\mathrm{mA} \, \int_0^\infty e^{-t/10} \, dt = 1 \,\mathrm{C}$$

$$N = \frac{1 \,\mathrm{C}}{1.6 \times 10^{-19} \,\mathrm{C/atom}} = 6.25 \times 10^8 \,\mathrm{atoms}$$


$$N = \frac{6.25 \times 10^{18}}{6.03 \times 10^{23}} = 1.036 \times 10^{-5} \,\mathrm{mol}$$


$$142.3 \times 1.036 \times 10^{-5} = 1.47 \times 10^{-3} \,\mathrm{g}$$


$$[\mathrm{Ag}^+][\mathrm{Cl}^-] = 1.56 \times 10^{-10}$$


$$[\mathrm{Ag}^+] = 1.73 \times 10^{-10} \,\mathrm{mol/liter}$$

$$1.73 \times 10^{-10} \times 142.3 = 2.46 \times 10^{-8} \,\mathrm{g}$$













## **Example**


We want to develop an electrical model for a specific biopotential electrode studies in the laboratory. The electrode is characterized by placing it in a physiological saline bath in the laboratory, along with an Ag/AgCl electrode having a much greater surface area and a known half-cell potential of 0.233 V. The dc voltage between the two electrodes is measured with a very-high-impedance voltmeter and found to be 0.572 V with the test electrode negative The magnitude of the impedance between two electrodes is measured as a function of frequency at very low currents; it is found to be that given in Figure in slide 12. From these data, determine a circuit model for the electrode.
























































Figure shows equivalent circuit of a biopotential electrode. A pair of these electrodes are tested in a beaker of physiological saline solution. The test consists of measuring the magnitude of the impedance between the electrodes as a function of frequency via low-level sinusoidal excitation so that the impedances are not affected by the current crossing the electrode-electrolyte interface. The impedance of the saline solution is small enough to be neglected. Sketch a Bode plot (log of impedance magnitude versus log of frequency) of the impedance between the electrodes over a frequency range of 1 to 100,000 Hz.

## Assume Figure in previous slide models both electrodes of the pair. The low corner frequency is Fc = $1/(2\pi RC) = 1/(2\pi \cdot 20 \text{ k}\Omega \cdot 100 \text{ n}F)$ = 80 Hz. The high corner frequency is Fc = $1/(2\pi \cdot 20 \text{ k}\Omega \cdot 100 \text{ n}F)$ = 5380 Hz. The slope between the two corner frequencies is –1 on a log-log plot.

## **Example**

A pair of biopotential electrodes are implanted in an animal to measure the electrocardiogram for a radiotelemetry system. One must know the equivalent circuit for these electrodes in order to design the optimal input circuit for the telemetry system. Measurements made on the pair of electrodes have shown that the polarization capacitance for the pair is 200 nF and that the half-cell potential for each electrode is 223 mV. The magnitude of the impedance between the two electrodes was measured via sinusoidal excitation at several different frequencies. The results of this information, draw an equivalent circuit for the electrode pair. State what each component in your circuit represents physically, and give its value.

Frequency	Impedance (Magnitude) ( $\Omega$ )
5 Hz	20,000
10 Hz	19,998
	*
	€
40 kHz	602
50 kHz	600
100 kHz	600

